



Quotation date: _____
Consultation was carried out by
Staff at the medical supplies shop: _____

Order date: _____
according to GTC PRO ACTIV Reha-Technik GmbH

Medical suppliers / company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Different invoice address (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____

Different shipping address (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____

Customer name / commission:

Surname, forename: _____
Street: _____
Postcode / city: _____
Phone: _____
E-Mail: _____
D.O.B: _____
Paying organisation: _____
Indication: _____

SPIKE adaptive bike

approved up to 75 kg trailer load and 5 kg payload

- docking plate for quick coupling to the "Quick-Snap" adapter
- drive unit with steering fork
- A-shaped handles with integrated switching and brake fittings
- chain guard in black on both sides
- without gear shift
- bottom bracket support can be adjusted in height and angle
- parking stand
- Sigma lighting set with rechargeable battery charger, front and rear lights

A9114000000 € 2,600.00

Basic equipment:

- Freewheel hub with singlespeed kit, 15 teeth
- Hydraulic disc brake with Ø 180 mm disc
- Mech. Shimano disc brake, brake handle left

Body size

Body height: _____ cm
Upper body height: _____ cm
Weight: _____ kg

This information are required for the processing.

Fork color

Information on current colour designs and example images at: www.proactiv-gmbh.de/farben.php

White RAL 9016

gloss A9000100170-023

matt A9000100170-123 € 42.10

Brilliant silver

gloss A9000100170-011

Anthracite RAL 7016

gloss A9000100170-001

matt A9000100170-101 € 42.10

Black RAL 9005

gloss A9000100170-003

matt A9000100170-019 € 42.10

Deep black effect silver

gloss A9000100170-007

Pure orange RAL 2004

gloss A9000100170-022

Fire red

gloss A9000100170-018

matt A9000100170-118 € 42.10

Pink magenta RAL 4010

gloss A9000100170-037

matt A9000100170-137 € 42.10

Sky blue

gloss A9000100170-017

Pastel blue RAL 5024

gloss A9000100170-036

matt A9000100170-136 € 42.10

Navy blue

gloss A9000100170-045

matt A9000100170-145 € 42.10

Electric green

gloss A9000100170-015

matt A9000100170-115 € 42.10

Lemon green

gloss A9000100170-005

matt A9000100170-105 € 42.10

Violet

gloss A9000100170-044

Yellow RAL 1023

gloss A9000100170-046

matt A9000100170-146 € 42.10

Red RAL 3001

gloss A9000100170-047

matt A9000100170-147 € 42.10

Lagoon blue

gloss A9000100170-052

matt A9000100170-152 € 42.10

Athletic grey

gloss A9000100170-050

matt A9000100170-150 € 42.10

Muddy grey

gloss A9000100170-051

matt A9000100170-151 € 42.10

Equipment

Shimano Deore 1x11 Rapid Fire 11-speed derailleur

A9114200032 € 391.90

mechanical dérailleur gear with operation on the handle
11-speed cassette 11-42 teeth

Shimano „Nexus“ 7-speed hub gear Twist shifter „Grip Shift“

A9114000007 € 527.00

with V-Brake rim brake, brake handle left



Shimano „Nexus“ 8-speed hub gear Twist shifter „Grip Shift“

A9114000008 € 697.80

with V-Brake rim brake, brake handle left



Mountaindrive 2.5 : 1 hill reduction in conjunction with hub gears

A9100200003 € 774.20


Planetary gear set for doubling the gears.
All gears (hub gears and hill reduction) can be shifted when stationary!
In conjunction with handlebar width 40 or 45 cm



Dimensions

Bottom bracket shaft with Standard crank

- Grip width 36 cm
A9114500100-036
- Grip width 40 cm
A9110500100-040
- Grip width 45 cm
A9110500100-045



Standard pedal crank

- Length 135 mm
A9110600100-135
- Length 155 mm
A9110600100-155
- Length 175 mm
A9110600100-175

Tyres

20x2" Schwalbe „Big Apple“

- A9114105125
(ETRTO 406-50)
Balloon tyres



20x1.75" Schwalbe „Marathon“

- A9110105121
(ETRTO 406-47)
Touring tyres



Puncture protection band for pneumatic tyres

- A9000206100-K **€ 54.20**

Accessories depending on wheelchair

Extended bike plate

- For SPEEDY A1 / A2
A9110100100-002 **€ 359.20**
- For BUDDY
A9110100100-003 **€ 359.20**
- For TRAVELER
A9110100100-001 **€ 359.20**
4° camber integrated
Track independently adjustable



Wheelbase extension via "locking mechanism with pull operation"

- For SPEEDY 4all /
4teen / 4you /
TRAVELER 4you **€ 548.30**
A9110102001
removable, with continuous axle tube
4° / 6° camber integrated
Track independently adjustable



Frame adapter

SPEEDY rigid frame adapter with "Quick-Snap-Adaption"

- A9110700141 **€ 502.40**
with quick-release clamp and centering pin holder
Colour: black



TRAVELER folding frame adapter with "Quick-Snap-Adaption"

- A9110701300 **€ 537.20**
with quick-release clamp and centering pin holder on the folding mechanism
Colour: black



BUDDY Starrahmenadapter mit „Quick-Snap-Adaption“

- A9110701400 **€ 502.40**
with quick-release clamp and centering pin holder
Colour: black



Wheelchair front adapter with universal clamp and "Quick-Snap-Adaption"

A9110701500 € 572.30
 adapter attachment to the wheelchair stem
 Operation on the right-hand side

Mounting situation for PA products:
 Straight front end from SB 31 cm to 54 cm
 V-shaped front end from SB 30 cm to 50 cm

PLEASE NOTE:
 If the exact tube diameter is not specified, the standard adapter for a tube diameter of 23-27 mm will be delivered.



Please select below:

A9110701525
 Tube diameter 23-27 mm

A9110701530
 Tube diameter 28-32 mm

A9110701535
 Tube diameter 33-37 mm

Clamp sleeve [70mm] included in standard
 For close distance between footrest and drive wheel

Clamp sleeve [110mm] € 44.50
 A9199000117
 For medium distance between footrest and drive wheel

Clamp sleeve [170mm] € 53.90
 A9199000154
 For a wide distance between footrest and drive wheel

Clamp lever with 3-way star handle for actuating the frame clamping

A9199701507 € 103.00
 mounted on the right and left for easy terminal operation with limited manual function
 Only in conjunction with front adapter



Accessories

Plug-in mudguard

out of plastic
 incl. required assembly components

A9122502100 € 106.20

In conjunction with luggage rack

A9122502101 € 56.40



Mudguard aluminum

incl. required assembly components

A9122500005 € 149.00

In conjunction with luggage rack

A9122500006 € 109.10

SIGMA BC 8.0 ATS bike computer

A9110901000 € 68.10
 wireless, with holder



Chain guard disc in frame colour

A9114500005 € 83.30

Aluminum luggage rack „cargo systems“


A9120900000 € 203.90

Colour: Black
 Max. load weight 10 kg
 max. 15 kg in conjunction with luggage bags with lower load height



Maneuvering castors for parking stand

A9110900017 € 118.50
 Maneuvering castors mounted on the parking stand via inlet



Options for users with reduced functions

Tetra hand support, pivotable, adjustable grip distance

A9110601101 € 164.50



Tetra shell handles

A9110601100 € 167.40



Quick release clip lever extension for rigid frame adapter

A9110700004 € 100.90

to make the opening and closing of the quick release clips easier.
only in conjunction with rigid frame adapter



Eccentric connection lever extension for Quick Snap connection


simplifies opening or closing of the Quick Snap Adapter. Tetra version

A9110700005 € 71.70

straight version
only in conjunction with rigid frame adapter

A9110700006 € 71.70

curved version
only in conjunction with front adapter



Adaptation help for centering pin

A9110700002 € 53.40

to simplify the insertion of the central adapter.



Chest strap for upper body fixation

size S (105 cm) A9000800016 € 231.80

size M (125 cm) A9000800013 € 231.80

size L (145 cm) A9000800014 € 231.80

size XL (165 cm) A9000800017 € 231.80

consisting of elastic webbing material in combination with neoprene with Velcro and operating tab.
Colour: black



Services

Delivery and instruction of the handbike / traction device by the specialist dealer incl. assembly of the frame adapter

on a PRO ACTIV wheelchair

A9990000008

prod. no: _____

on a third-party product

A9990000005

Manufacturer: _____

Model: _____

Frame tube ø ____ mm

Notes and remarks:

Location, Date

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

- 5 - 2024-06-01 / The prices are plus statutory value-added tax and shipping costs

= Included in the basic price = Position selected

Order sheet Spike adaptive bike Rev. 06 UK

Signature End customer/dealer