



TRAVELER 4all Ergo

SIT ERGONOMICALLY,
TRAVEL EFFORTLESSLY.

The TRAVELER 4all Ergo active folding wheelchair features fine frame geometry with an ergonomic seat fold. You can choose 9 frame parameters to tailor it to your exact needs. With a compact pack size of 17 x 35 x 58 cm, it even fits in a sports bag. Thanks to the minimised frame geometry, the proven folding mechanism, the easy-to-fold backrest and the ability to remove assemblies without tools, it is also ideal for independent loading.




#Reliability

Thanks to its excellent stability, you can combine the TRAVELER 4all Ergo with adaptive bikes, traction devices and steering attachments as well as other adaptations.



#Flexibility

Many settings, such as tipping point, seat height, seat angle or back angle, can also be adjusted to your needs afterwards.



#Independence

With its small pack size, the folding backrest and tool-free detachable assemblies, the TRAVELER 4all Ergo is easy to transport.



#MoreMobility

By simply folding it, you can reduce the seat width (sitting in a wheelchair) by 15% or 5-10 cm to be able to pass through narrow passages effortlessly.





PERFECTLY ADAPTED TO YOUR LIFE.

Your **TRAVELER 4all Ergo** is manufactured entirely according to your wishes, from the frame dimensions to the equipment. You can choose from options such as footrests, brakes, push handles, back and seat systems, clothing guards armrests, and drive and caster wheels to tailor your TRAVELER 4all Ergo perfectly to your needs.

UNIQUE, SELECTABLE OPTIONS FROM PRO ACTIV:



The 7-way angle-adjustable and foldable backrest offers maximum flexibility for the ideal sitting position.



The footrest, which can be folded towards the back with a spring locking mechanism, is easy to handle and flexible in everyday life.



Elegant integral parking brake with comfortable and safe operation.



The aluminium rim in an elegant 3-spoke design combines maximum stability with an ultra-light weight – for the highest demands without compromise.

TECHNICAL SPECIFICATIONS

Seat width:	33 - 52 cm (in 1 cm steps)
Seat depth:	36 - 48 cm (in 1 cm steps)
Back height:	20 - 48 cm (in 1 cm steps)
Lower leg length:	37 - 51 cm (in 1 cm steps)
Rear seat height:	39 - 45 cm (in 1 cm steps)
Front seat height:	44 - 52 cm (in 1 cm steps)
Tipping point:	3 - 11 cm (in 1 cm steps)
Total product weight:	Starting at 11.2 kg
Max. standard load weight:	Up to 120 kg

Data subject to design changes.

FOLLOW US ON SOCIAL MEDIA



PROACTIV®

PRO ACTIV Reha-Technik GmbH

Im Hofstätt 11
72359 Dotternhausen
info@proactiv-gmbh.de

www.proactiv-gmbh.com

Rev. 09/2024

Subject to technical modifications, printing errors & colour deviations.