

Im Hofstätt 11
D-72359 Dotternhausen
Tel +49 (0)7427 9480-0
Fax +49 (0)7427 9480-7025
e-mail info@proactiv-gmbh.de
www.proactiv-gmbh.de



FREAK compactbike

Maximum user weight limit 100 kg

- Underframe chassis with reinforcement strut
- chassis length 102 cm, SB 35 cm
- Drive unit can be uncoupled via excenter lock
- black drive unit with 20" drive wheel
- A-shaped handles with integrated gearshift controls
- Shimano V rim brake
- Backrest with continuous angle adjustment
- wheel camber 6°, rear wheel size 20"
- safety flag
- Seat upholstery Body Contour

A9113000013

€ 4,640.00

Quotation date: _____

Consultation was carried out by
Staff at the medical supplies shop: _____

Order date: _____

according to GTC PRO ACTIV Reha-Technik GmbH

Medical suppliers / company address:

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Phone: _____

Fax: _____

E-Mail: _____

Contact: _____

Different invoice address (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Different shipping address (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Customer name / commission:

Surname, forename: _____

Street: _____

Postcode / city: _____

Phone: _____

E-Mail: _____

D.O.B: _____

Paying organisation: _____

Indication: _____

Body size

Body height: _____ cm

Weight: _____ kg

Frame colour - one-colour design

Information on current colour designs and example images at: www.proactiv-gmbh.de/farben.php

White RAL 9016

- gloss A9000100170-023
 matt A9000100170-123 € 42.10

Brilliant silver

- gloss A9000100170-011

Anthracite RAL 7016

- gloss A9000100170-001
 matt A9000100170-101 € 42.10

Black RAL 9005

- gloss A9000100170-003
 matt A9000100170-019 € 42.10

Deep black effect silver

- gloss A9000100170-007

Pure orange RAL 2004

- gloss A9000100170-022

Fire red

- gloss A9000100170-018
 matt A9000100170-118 € 42.10

Pink magenta RAL 4010

- gloss A9000100170-037
 matt A9000100170-137 € 42.10

Sky blue

- gloss A9000100170-017

Pastel blue RAL 5024

- gloss A9000100170-036
 matt A9000100170-136 € 42.10

Navy blue

- gloss A9000100170-045
 matt A9000100170-145 € 42.10

Electric green

- gloss A9000100170-015
 matt A9000100170-115 € 42.10

Lemon green

- gloss A9000100170-005
 matt A9000100170-105 € 42.10

Violet

- gloss A9000100170-044

Yellow RAL 1023

- gloss A9000100170-046
 matt A9000100170-146 € 42.10

Red RAL 3001

- gloss A9000100170-047
 matt A9000100170-147 € 42.10

Lagoon blue

- gloss A9000100170-052
 matt A9000100170-152 € 42.10

Athletic grey

- gloss A9000100170-050
 matt A9000100170-150 € 42.10

Muddy grey

- gloss A9000100170-051
 matt A9000100170-151 € 42.10

Neodrives

e-assistant version with neodrives

- A9113200000 € 3,155.40

- Neodrives Z10 motor
 - torque 12/40 Nm, power 250 W
 - colour display sMMi Z15 Connect
 - push assist for riding without pedal assistance up to 6 km/h
 - Neodrives RR G2 battery with 14.25 Ah / 513 Wh
 - range approx. 40 - 70 km
 Not in conjunction with hub gears!



Hub gear

Not in conjunction with e-assistant version

Shimano „Nexus“ 8- speed hub gear

- A9113000015
 „Grip Shift“ twist shift
 Gear ratio: 32 teeth sprocket,
 cassette 20 teeth
 with back-pedal brake
 Shimano V-Brake brake handle left

Shimano „Alfine“ 11- speed hub gear

- A9113200026 € 809.40
 „Rapid Fire“ gear lever
 Gear ratio: 32 teeth sprocket,
 cassette 20 teeth
 with hydraulic disc brake (Magura MT2) brake
 handle right and Shimano V-Brake brake
 handle left



Rohloff „Speedhub 500/14“ 14- speed hub gear

- A9113000016 € 1,712.60
 „Grip Shift“ twist shift
 Gear ratio: 32 teeth sprocket,
 cassette 16 teeth
 with hydraulic disc brake (Magura MT2) brake
 handle right and Shimano V-Brake brake
 handle left

derailleur gear

Shimano Deore 1x11 Rapid Fire 11-speed derailleur

A9113200034 € 502.40

mechanical derailleur gear with operation on the handle
11-speed cassette 11-42 teeth
with hydraulic disc brake (Magura MT2) brake handle right and Shimano V-Brake brake handle left
Top chainring with chain guards

Selection of the chainring:

- 32 teeth
 38 teeth

Shimano Deore XT 1x11 DI2 with push button operation 11-speed derailleur

A9113200035 € 1,747.80

electronic derailleur gear with control on the handle. 11-speed cassette 11-42 teeth
with hydraulic disc brake (Magura MT2) brake handle right and Shimano V-Brake brake handle left
Top chainring with chain guards

Selection of the chainring:

- 32 teeth
 38 teeth



Shimano Deore XT 1x11 DI2 with chin control 11-speed derailleur

A9113200033 € 1,900.30

electronic derailleur gear with chin control (tetra equipment).
11-speed cassette 11-42 teeth
with hydraulic disc brake (Magura MT2) brake handle right and Shimano V-Brake brake handle left
Top chainring with chain guards

Selection of the chainring:

- 32 teeth
 38 teeth



Gearshift accessories

Mountaindrive 2.5 : 1 Hill reduction

A9100200003 € 774.20

Planetary gear set for doubling of the number of gears permits switching to the gear reduction for hills when stationary.
As standard with 38 tooth chainring.
Only in conjunction with grip width 40 or 45 cm



Brake system

Powerful braking system with brake disk Ø 203 mm

A9110300105 € 80.70

Recommended from 80 kg bodyweight



Dimensions

Crank type and grip width

Bottom bracket shaft with Standard crank

- Grip width 40 cm
A9110500100-040
- Grip width 45 cm
A9110500100-045

Bottom bracket shaft with Tetra crank (flared)

with extended crank arms for maximum upper body space

- Grip width 48 cm (176 mm)
A9110500100-048
- Grip width 53 cm (226 mm)
A9110500100-053

Standard pedal crank

- Length 155 mm
A9110600100-155
- Length 175 mm
A9110600100-175
- Length 195 mm
A9110600100-195



Tetra pedal crank


- Length 155 mm € 132.10
A9100600105-155
- Length 175 mm € 132.10
A9100600105-175
- individual length
140 mm to 195 mm € 251.20
A9100600100



Fixed bottom bracket support

*(Limited availability - until stocks last.
Alternative: adjustable bottom bracket support)*

- S (X/Y 67/153 mm)
A9100501100-003



Adjustable bottom bracket support

adjustable in angle via clamped pivot point at the fork bridge and in length

- S (180 mm) € 225.30
A9100500001-002
- M (240 mm) € 225.30
A9100500001-003

Please note: If not specified, the adjustable bottom bracket support will be delivered in M.



Legrest

Divided legrest

- A9100100008
adaptation on the fork
for short legs
Only possible up to 165 cm body height



Foot stop for legrest

- A9100100012 € 88.50



Upper thigh fixation

- A9100100013 € 208.80



Back equipment

Backrest with filler panel and back cushion

A9113800002
Cover: Black, breathable, open surface

Back upholstery Body Contour

A9113800010 € 86.80
with back cushion 3 cm
Cover: Black, breathable, open surface



Seat cushion

Seat cushion „Standard“

A9113000190 € 108.70
for chassis 102 cm
SB: 35 cm
height: 4 cm
Cover: Black, breathable, open surface



Tyres

20x1.1" Schwalbe „Durano“

A9113105001
(ETRTO: 406-28)
racing bike tyres



20x1.75" Schwalbe „Marathon Plus“

A9113105120 € 69.00
(ETRTO: 406-47)
unpuncturable
with anti-puncture belt Smart Guard



Accessories

Plug-in mudguard

A9100502100 € 56.40



Leg guard

A9100100231 € 166.00
protection for MTB tires
bolted to the fork leg
Colour: black matt

Collision guard clamped on axle tube

A9100800100 € 211.70
with mounting clamp 40 mm



Neck rest

A9100800150 € 253.40
height-adjustable
angle-adjustable



Pelvic fixation belt

A9100800200 € 115.20



SIGMA BC 8.0 ATS bike computer

A9110901000 € 68.10
wireless, with bracket
Not in conjunction with e-assistant version



SIGMA BC 14.0 WL CAD bike computer

A9110900018 € 166.50

with cadence set, wireless
Not in conjunction with e-assistant version




Bracket ø 22 mm x 42 mm for mounting accessories (tacho, lighting, etc.)

A9110903271 pcs. á € 22.50

Colour: natural anodized

_____ piece



SIGMA lighting set

A9120900019 € 120.30

consisting of front light and rear lights, chargeable via Micro USB socket, incl. Micro USB charging cable. Lighting set compliant with road traffic regulations



Options for users with reduced functions

Tetra hand support, pivotable, adjustable grip distance

A9110601101 € 164.50



Tetra shell handles

A9110601100 € 167.40



Tetra backrest swivel and adjustable grip width

A9110601102 € 186.00

very stiff external support as round tube ø 20 mm with foam rubber coating



Services

Delivery and instruction of the handbike / traction device by the specialist dealer incl. assembly of the frame adapter

on a PRO ACTIV wheelchair

A9990000008

prod. no: _____

on a third-party product

A9990000005

Manufacturer: _____

Model: _____

Frame tube ø _____ mm

Notes and remarks:

Location, Date

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

- 6 - 2024-06-01 / The prices are plus statutory value-added tax and shipping costs

= Included in the basic price = Position selected