



PRO ACTIV Reha-Technik GmbH
 Im Hofstätt 11
 D-72359 Dotternhausen
 Phone +49 (0)7427 9480-0
 Fax +49 (0)7427 9480-7025
 E-mail: info@proactiv-gmbh.de
 www.proactiv-gmbh.com



☐ Offer date: _____

☐ Order date: _____
 according to GTC PRO ACTIV Reha-Technik GmbH

Consultant: _____

*** Required fields, please fill in completely!**

Order address / customer:

Cust. no.:* _____

Last name:* _____

First name:* _____

Street:* _____

Postcode / city:* _____

Phone:* _____

Fax: _____

E-mail:* _____

Collection: ☐ Yes ☐ No

Delivery address

Name: _____

Street: _____

Postcode / city: _____

Phone: _____

Fax: _____

E-Mail: _____

Date of birth: _____

Note: _____

Biometric data:

Body height: _____

Weight: _____

Inside leg: _____

PED-E

€ 4,142.02

- Aluminium frame with CNC milled/welded construction
- Aluminium fork with milled fork bridges,
Colour: black
- Hydraulic disc brakes 180 mm / 160 mm
- Neodrives Z15 drive system with TFT display
- Rechargeable battery integrated in the frame
37 V; 16 Ah; 587 Wh
- Shimano Zee dérailleur 1x10
- Pedal crank (L = 80 mm)
- Aluminium rims
- Saddle support Ø 27.2 mm
- Side stand
- Bell
- Lighting set
- Cargo load: 75 kg

Wheel size

16"

☐ A9525000160

Schwalbe Big Apple 16x2"

☐ A9525200005
(50-305)
Performance Line
RaceGuard E25
Puncture protection level 4
Colour: black



18"

☐ A9525000180

Schwalbe Big Apple 18x2"

☐ A9525200010
(50-355)
Performance Line
RaceGuard E25
Puncture protection level 4
Colour: black



20"

☐ A9525000200


Schwalbe "Jumpin' Jack" 20x2.1"

☐ A9525105100
(54-406)
Performance Line
Puncture protection level 3.5
Colour: black



Schwalbe "Smart Sam" 20x2.35"

☐ A9525105101 € 24.79
(60-406)
Performance Line
Puncture protection level 3.5
Colour: black



Schwalbe "Big Ben Plus" 20x2.15"


☐ A9525105102 € 24.79
(55-406)
Performance Line
DD, GreenGuard E50
Puncture protection level 6
Colour: black



Bottom bracket position

Please specify the desired bottom bracket position:

If nothing is specified, the standard bottom bracket position is 1.



Saddle

Velo Fit Traveller saddle

☐ A9500800010

Size M

Width: 15 cm; length: 27 cm

☐ A9500800020

Size L

Width: 19 cm; length: 26 cm

Colour: black



Definition of the saddle support length

Distance from the upper edge of the rechargeable battery to the upper edge of the saddle:

_____ cm

If nothing is specified, the saddle support will be delivered with extended length (40 cm) and can be adjusted by the customers themselves.



Front end

Front end 7°

☐ Length: 90 mm
A9500400500-090

☐ Length: 110 mm
A9500400500-110

Clamping: Oversize 31.8 mm

Colour: black



Ergotec Swell RE front end

☐ A9500400550 **€ 39.58**

Length: 100 mm

Clamping: Oversize 31.8 mm
angle-adjustable

Colour: black



Handlebar

Sport Ergotec handlebar

☐ A9500400610

Width: 720 mm

Height: 15 mm

Grip bend: 9°

Oversize 31.8 mm

Colour: black

The handlebar width can be individually adjusted in the range from 720 to 615 mm:
_____ mm



Comfort Ergotec handlebar

☐ A9500400620

Width: 680 mm

Height: 30 mm

Grip bend: 14°

Oversize 31.8 mm

Colour: black

The handlebar width can be individually adjusted in the range from 680 to 615 mm:
_____ mm



Grips

Sport grips black / grey

☐ A9500400700

Shock-absorbing gel grip surface,
anti-twist clamping on the handlebar



Gel Ergo Wing grip cover

☐ A9500400710 **€ 21.26**

Ergonomic gel grip surface,
anti-twist clamping on the handlebar



Gel Ergo Wing Comfort grip cover

☐ A9500400720 **€ 26.81**

Large, ergonomically shaped wrist rest,
anti-twist clamping on the handlebar



Pedals

Offroad pedals

☒ A9500500300
CR-MO axle,
with reflector (K-mark)
Platform with jagged profile



Touring pedals

☒ A9500500310
CR-MO axle,
with reflector (K-mark),
rubberised support



Ergo Grip pedals

☐ A9500500320
CR-MO axle,
with reflector (K-mark),
sandpaper support

€ 26.13



Accessories

Mudguard 53 mm, Colour-coated on front and back

☐ A9525700053 € 191.43
not for 16"
Colour selection:
see colour design at the back



Luggage rack

☐ A9525700100 € 416.47
Aluminium
Incl. rear light with battery



Colour design – frame and attachment parts

For the following components, a colour can be chosen from the information sheet “Pedelec Colour Design”:

RAL powder paint or “Dormant” powder paint or “Liquid coating” custom paint

Frame colour:

☒ RAL powder paint and “Dormant” powder paint:

Please enter the last 3 digits of the colour:

A9500100170 - in words (colour): _____

☐ “Liquid coating” custom paint:

Please enter the last 2 digits of the custom colour:

A9500100170 - 5 in words (colour): _____

€ 161.68

Mudguard colour:

☒ RAL powder paint or “Dormant” powder paint:

A9500700170

Please enter the last 3 digits of the colour:

A9500100170 - in words (colour): _____

☐ “Liquid coating” custom paint:

A9500700172

Please enter the last 2 digits of the custom colour:

A9500100170 - 5 in words (colour): _____

€ 111.85

Notes and remarks:

I hereby confirm the information provided above and agree to receiving a written offer without obligation.

Location, date End customer signature (offer)

Location, date Consultant signature (offer)

Location, date End customer signature (order)

Location, date Consultant signature (order)

Info sheet: Colour Designs Pedelecs

Colour Designs - Powder Coatings

Information on current colour designs and sample images at: www.proactiv-gmbh.com/colours/

Powder Coatings

white RAL 9016

☒ gloss A9500100170-023

anthracite RAL 7016

☒ gloss A9500100170-001

☒ matt A9500100170-101

black RAL 9005

☒ gloss A9500100170-003

☒ matt A9500100170-019

pure orange RAL 2004

☒ gloss A9500100170-022

☒ matt A9500100170-122

yellow RAL 1023

☒ gloss A9500100170-046

☒ matt A9500100170-146

red RAL 3001

☒ gloss A9500100170-047

☒ matt A9500100170-147

„Dormant“ Powder Coatings

sky blue

☒ gloss A9500100170-017

electric green

☒ gloss A9500100170-015

violet

☒ gloss A9500100170-044

fire red

☒ gloss A9500100170-018

Powder Coatings “Trend colours”

athletic grey

☒ gloss A9500100170-303 € 30.59

☒ matt A9500100170-304 € 30.59

muddy grey

☒ gloss A9500100170-305 € 30.59

☒ matt A9500100170-306 € 30.59

lagoon blue

☒ gloss A9500100170-301 € 30.59

☒ matt A9500100170-302 € 30.59

Colour Designs - Custom Colours ‘Liquid Coating’

€ 161.68

storm (ice blue)

☐ gloss A9500100170-500

☐ matt A9500100170-501

coffee

☐ gloss A9500100170-520

☐ matt A9500100170-521

prairie (rose gold)

☐ gloss A9500100170-504

☐ matt A9500100170-505

blue ocean

☐ gloss A9500100170-522

☐ matt A9500100170-523

mint

☐ gloss A9500100170-528

☐ matt A9500100170-529