



STABILITY, THAT CAN BE FOLDED TOGETHER.

As great as the performance of the TRAVELER is in daily use, it is just as small when it is to be transported. You can load the wheelchair into the car and unload it again yourself, offering you a high level of independence. With its extremely torsion-resistant and reliable folding mechanics, the TRAVELER combines the handling characteristics of a rigid frame wheelchair with the advantages of a folding wheelchair.

TRAVELER



#Reliability

Thanks to its excellent stability, you can combine the TRAVELER with adaptive bikes, towing and attachments as well as other adaptations.



#Flexibility

Many settings, such as tipping point, seat height, seat angle or back angle, can also be adjusted to your needs afterwards.



#Independence

With its small pack size, folding backrest and tool-free detachable assemblies, the TRAVELER is easy to transport.



#MoreMobility

By simply folding it, you can reduce the seat width (sitting in a wheelchair) by 15% or 5-10 cm to be able to pass through narrow passages effortlessly.



FLEXIBLE, STURDY AND ALWAYS GOOD FOR A TRIP.

Your **TRAVELER** is individually manufactured according to your wishes, from the frame dimensions to the equipment. You can choose from options such as footrests, brakes, push handles, seat and back systems, armrests, and drive and caster wheels to tailor the TRAVELER perfectly to your needs.


UNIQUE, SELECTABLE OPTIONS FROM PRO ACTIV:



The 7-way angle-adjustable and foldable backrest offers maximum flexibility for the ideal sitting position.



The footrest, which can be folded towards the back with a spring locking mechanism, is easy to handle and flexible in everyday life.



Manoeuvrable one-hand drive with double handrim for users with hemiplegia or restricted mobility on one side of the body.



With the Heavy Duty equipment, your TRAVELER is designed for high loads and is suitable for user weights up to 180 kg.



TECHNICAL SPECIFICATIONS

Seat width:	33 - 52 cm (in 1 cm steps)
Seat depth:	36 - 48 cm (in 1 cm steps)
Back height:	20 - 48 cm (in 1 cm steps)
Lower leg length:	37 - 51 cm (in 1 cm steps)
Rear seat height:	37 - 43 cm (in 1 cm steps)
Front seat height:	44 - 52 cm (in 1 cm steps)
Tipping point:	4 - 13 cm (in 1 cm steps)
Total product weight:	Starting at 11.5 kg
Max. standard load weight:	Up to 120 kg
Max. load weight for Heavy Duty:	Up to 180 kg appropriate individual dimensions upon request

Data subject to design changes.

FOLLOW US ON SOCIAL MEDIA



PROACTIV[®]

PRO ACTIV Reha-Technik GmbH

Im Hofstätt 11

72359 Dotternhausen

info@proactiv-gmbh.de

www.proactiv-gmbh.com

Rev. 09/2024

Subject to technical modifications, printing errors & colour deviations.