



Quotation date: _____
Consultation was carried out by
Staff at the medical supplies shop: _____

Order date: _____
according to GTC PRO ACTIV Reha-Technik GmbH

Medical suppliers / company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Different invoice address (see below)
Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____

Different shipping address (see below)
Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____

Customer name / commission:
Surname, forename: _____
Street: _____
Postcode / city: _____
Phone: _____
E-Mail: _____
D.O.B: _____
Paying organisation: _____
Indication: _____

HUSK-E adaptive bike

Comfortable version with Neodrives system up to 120 kg trailer load

- For drivers with low muscle strength and stronger impairments as well as for all those who need a comfortable drive with high power assistance
- For uneven terrain without long climbs
- with Pinion or derailleur gear
- Disc brake Ø 180 mm with brake handle or back-pedal disc brake possible

Standard equipment

- Neodrives Hub motor Z15 | 250 W | 40 Nm
- Docking plate for quick attachment to the "Quick-Snap-Adapter"
- adjustable and foldable bottom bracket support to minimize the packing size
- Hydraulic disc brake
- Mechanical parking disc brake
- Parking stand with longer operating lever
- Lighting set
- Mudguard
- 20" drive wheel
- Colour display sMMi Z15 connect

A9117100040-020 € 7,390.00

Body size

Body height: _____ cm
Upper body height: _____ cm
Weight: _____ kg

This information are required for the processing.

Transmission

Chainring on bottom bracket with 28 teeth

A9117200050
For a higher top speed

Chainring on bottom bracket with 22 teeth

A9117200051
For a very light response

Derailleur gears

Shimano Deore 1x11 Rapid Fire 11-speed derailleur gears

A9117200029
mechanical derailleur gear with control on the handle.
11-speed cassette 11-42 teeth and a gear range of 381 %
Gearshift button mounted on the right handle.



Shimano Deore XT 1x11 Rapid Fire 11-speed derailleur gears

A9117200019 € 216.80
mechanical derailleur gear with control on the handle.
11-speed cassette 11-42 teeth and a gear range of 381 %
Gearshift button mounted on the right handle.



Shimano Deore XT 1x11 Di2 with push button operation 11-speed derailleur gears

A9117200028 € 1,442.80
electronic derailleur with operation on the handle via button.
11-speed cassette 11-42 teeth and a gear range of 381 %

Mounting of the switching button:


Right side Left side

If not specified, the push-button will be delivered on the right side.



Shimano Deore XT 1x11 Di2 with chin control 11-speed derailleur gears

A9117200027 € 1,583.60
electronic derailleur gear with chin control (Tetra equipment).
11-speed cassette 11-42 teeth and a gear range of 381 %



Pinion gear shift

Pinion C1.6 gear shift

A9117200010 € 1,226.40
with 6 gears
transmission range of 295 %
Operation on the handle via a rotating switch.

Mounting: Right side Left side
If not specified, the mounting will be on the right side.



C1.6

Pinion C1.9 XR gear shift

A9117200011 € 1,368.90
with 9 gears
transmission range of 568 %
Operation on the handle via a rotating switch.

Mounting: Right side Left side
If not specified, the mounting will be on the right side.



C1.9 XR

Brake system

Double hydraulic disc brake system „Magura”

A9117300100 € 163.20

Two Magura 2-piston callipers
Parking brake via locking bar on the handle
Omitting the mechanical parking disk brake.

Mechanical back-pedal disc brake with crank release function

A9117300110 € 716.90

compatible with derailleur gears and freewheel hub. Incl. second mechanical disc brake designed as a parking brake.
The hydraulic disc brake is omitted.
In conjunction with:
- grip width 45 cm, 50 cm, 53 cm or 58 cm
- crank length 155 or 175 mm
up to 120 kg body weight



Driving brake - supplement for the mechanical parking brake

A9117300120 € 112.70

The parking disc brake is supplemented with a brake handle on the left-hand handle. This means that the mechanical disc brake can be used as a driving brake and parking brake.

Batteries and accessories

The specified ranges depend on the distance profile, user weight and speed.

Battery Neodrives DT 36V 16Ah 576Wh

Range approx. 45-75 km

A9117900140
Mounting on the left of the fork



Double Battery Neodrives DT 36V 16Ah 576Wh

Range approx. 90-150 km

A9117900145 € 1,425.20

Mounted on both sides of the fork, interchangeable connector mounted on the fork for easy switching between the batteries.
Bracket swivels to the right for removing the battery.
Only in conjunction with frame adapter with right-hand operation.



Colour mudguard

If not specified, the mudguard and chain guard will be delivered in anthracite RAL 7016 matt.

Mudguard, chain guard and motor cover Anthracite RAL 7016 matt

A9117500004

Mudguard, chain guard and motor cover Pure orange RAL 2004 gloss

A9117500006



Mudguard, chain guard and motor cover Lemon green gloss

A9117500007



Mudguard, chain guard and motor cover Customized colour selection

A9117500010 € 149.00

from colour palette one-colour design
powder colours wheelchair

A9000100170- _____

Colour: _____

Tyres

20x2.35" Schwalbe „Smart Sam“

A9117105101
(ETRTO 406-60)
profile tyre
puncture protection level 3.5



20x2.15" Schwalbe „Big Ben Plus“

A9117105102 € 27.90
(ETRTO 406-55)
profile tyre
puncture protection level 6



Dimensions

Crank type and grip width

Standard crank

- Grip width 40 cm
A9117500100-040
- Grip width 45 cm
A9117500100-045
- Grip width 50 cm
A9117500100-050

Tetra crank (flared)

with extended crank arms for maximum upper body space

- Grip width 48 cm
A9117500100-048
- Grip width 53 cm
A9117500100-053
- Grip width 58 cm
A9117500100-058

Bottom bracket support

maximum length from the pivot point in the fork to the bottom bracket shaft

- Length 520 mm (short)
using shortened holding tubes € 150.80
A9117500001
Recommendation: up to approx. 165 cm body height
- Length 570 mm (standard)
A9117500008
Recommendation: from approx. 165 cm to 185 cm body height
- Length 620 mm (long)
using a longer bottom bracket distance € 198.50
A9117500003
Recommendation: from approx. 185 cm body height

Standard pedal crank

- Length 155 mm
A9110600100-155
- Length 175 mm
A9110600100-175
- Length 195 mm
A9110600100-195



Tetra pedal crank

- Length 155 mm € 132.10
A9100600105-155
- Length 175 mm € 132.10
A9100600105-175
- individual length
140 mm to 195 mm € 251.20
A9100600100



Handles

Palm rest on the handle

- A9110601200 € 99.80
upholstered, bolted, angle-adjustable.



Tetra hand support, pivotable, adjustable grip distance

- A9110601101 € 164.50



Tetra backrest swivel and adjustable grip width

- A9110601102 € 186.00
very stiff external support as round tube
ø 20 mm with foam rubber coating






<p>Tetra ergonomic handle free-standing</p> <p><input type="checkbox"/> A9100600001 € 159.80</p> <p>Crank handles arranged vertically Can be combined with Tetra hand rest, pivotable</p> 	<p>Tetra shell handles</p> <p><input type="checkbox"/> A9110601100 € 167.40</p> 	
--	---	--

Accessories Bike

<p>Maneuvering castors for parking stand</p> <p><input type="checkbox"/> A9110900017 € 118.50</p> <p>Maneuvering castors mounted on the parking stand via inlet</p> 	<p>LED lighting Neodrives system</p> <p><input type="checkbox"/> A9117900205 € 147.30</p> <p>LED front light powered by the drive system battery Rear light as plug-in light with battery supply Lighting set from the basic equipment omitted</p> 	<p>Removable luggage rack</p> <p><input type="checkbox"/> A9117900210 € 315.60</p> <p>Colour: Black Max. load weight 10 kg</p> <p>Chain tensioner</p> <p><input type="checkbox"/> A9117200034 € 45.20</p> <p>Additionally required for the luggage rack if a hub gear or Pinion gear system is selected.</p> 
---	---	---

<p>Operating satellite on the handle</p> <p>Mounting: Right side <input type="checkbox"/> Left side <input type="checkbox"/></p> <p>If not specified, the mounting is on the left side.</p> <p>with existing support line (laid along a brake line)</p> <p><input type="checkbox"/> A9116200000 € 98.20</p> <p>without existing support line (an additional blind line is required)</p> <p><input type="checkbox"/> A9116200001 € 139.20</p>	
---	--

Options for users with reduced functions

<p>Quick release clip lever extension for rigid frame adapter</p> <p><input type="checkbox"/> A9110700004 € 100.90</p> <p>to make the opening and closing of the quick release clips easier. only in conjunction with rigid frame adapter</p> 	<p>Adaptation help for centering pin</p> <p><input type="checkbox"/> A9110700002 € 53.40</p> <p>to simplify the insertion of the central adapter.</p> 	<p>Lever extension excenter lock for Quick Snap connection</p> <p><input type="checkbox"/> A9110700006 € 71.70</p> <p>simplifies the opening and closing of the Quick Snap adapter Tetra version, curved</p> 
--	---	--

Chest strap for upper body fixation

- size S (105 cm) A9000800016 € 231.80
- size M (125 cm) A9000800013 € 231.80
- size L (145 cm) A9000800014 € 231.80
- size XL (165 cm) A9000800017 € 231.80

consisting of elastic webbing material in combination with neoprene with Velcro and operating tab.

Colour: black



Frame adapter

SPEEDY rigid frame adapter right-hand operation with "Quick-Snap-Adaption"

- A9110700142 € 502.40

with quick-release clamp and centering pin holder

Colour: black



TRAVELER folding frame adapter right-hand operation with "Quick-Snap-Adaption"

- A9110701302 € 537.20

with quick-release clamp and centering pin holder on the folding mechanism

Colour: black



BUDDY rigid frame adapter right-hand operation with "Quick-Snap-Adaption"

- A9110701401 € 502.40

with quick-release clamp and centering pin holder

Colour: black



Wheelchair front adapter right-hand operation with universal clamp and "Quick-Snap-Adaption"

- A9110701502 € 572.30

adapter attachment to the wheelchair stem
Operation on the right-hand side

Mounting situation for PA products:
Straight front end from SB 31 cm to 54 cm
V-shaped front end from SB 30 cm to 50 cm

PLEASE NOTE:

If the exact tube diameter is not specified, the standard adapter for a tube diameter of 23-27 mm will be delivered.



Please select below:

- A9110701525

Tube diameter 23-27 mm

- A9110701530

Tube diameter 28-32 mm

- A9110701535

Tube diameter 33-37 mm

- Clamp sleeve [70mm] included in standard

For close distance between footrest and drive wheel

- Clamp sleeve [110mm] € 44.50 A9199000117

For medium distance between footrest and drive wheel

- Clamp sleeve [170mm] € 53.90 A9199000154

For a wide distance between footrest and drive wheel

Clamp lever with 3-way star handle for actuating the frame clamping

A9199701507 € 103.00

mounted on the right and left for easy terminal operation with limited manual function
Only in conjunction with front adapter



Accessories depending on wheelchair


Extended bike plate

For SPEEDY A1 / A2
A9110100100-002 € 359.20

For BUDDY
A9110100100-003 € 359.20

For TRAVELER
A9110100100-001 € 359.20

4° camber integrated
Track independently adjustable



Wheelbase extension via „locking mechanism with pull operation”

For SPEEDY 4all /
4teen / 4you /
TRAVELER 4you € 548.30
A9110102001

removable, with continuous axle tube
4° / 6° camber integrated
Track independently adjustable



Wheelbase extension angle-adjustable and removable

For SPEEDY F2 / F4
A9015100005 € 711.20

with continuous axle tube,
height adjustment via toothed joint
adjustable in 6° angle increments
Wheel position displacement approx. 14 cm
4° / 6° camber integrated, toe adjustable
Mounting option must be checked by PRO
ACTIV.

Only possible in conjunction with base support
plate.
up to 120 kg trailer load.



Services

**Delivery and instruction of the
handbike / traction device by the
specialist dealer incl. assembly of the
frame adapter**

on a PRO ACTIV wheelchair

A9990000008

prod. no: _____

on a third-party product

A9990000005

Manufacturer: _____

Model: _____

Frame tube ø ____ mm

Notes and remarks:

Location, Date

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

- 7 - 2024-06-01 / The prices are plus statutory value-added tax and shipping costs

= Included in the basic price = Position selected

Order sheet HUSK-E adaptive bike Z15 Rev. 01 UK